

Safety + Warranty

MOBILE COMPANION 2

Mobile Companion 2: An on-the-go safety communication accessory for independent adults. **Included:** Aloe Care Mobile Companion, charging station, power cord, plug adapter, belt clip, and lanyard. **Required (not included):** Android or iOS mobile device.

The Mobile Companion and other accessory devices, including but not limited to the Aloe Care ecosystem (collectively, "Equipment"), are provided "as is," and Aloe Care Health, Inc. ("Aloe Care") hereby disclaims to the fullest extent permitted by law, all liability in connection with the equipment and the use thereof, including without limitation in connection with claims related to personal injury, personal property, or the infringement, misappropriation, or other violation of intellectual property rights. Aloe Care also hereby makes no guarantees related to the equipment and disclaims any warranties, whether written or oral or express or implied, including any warranty of quality, merchantability, or fitness for a particular use or purpose in connection with the equipment.

Fall detection does not detect 100% of falls. The equipment and medical alert services are not a substitute for 911. If you need help and can dial 911, you should do so.

IPX7: Do not use or store in environments prone to dust, moisture, or extreme temperatures. Device protected against water projected in powerful jets (12.5 mm (0.49 in)) against the enclosure from any direction and short periods of immersion in water—between 15 cm (5.9 inches) and 1 meter (3.3 feet) for up to 30 minutes.

Warning: Choking and strangulation hazard—the equipment contains small parts and cords that may present a choking or strangulation risk for children and pets. This equipment is not intended for children under 12 years of age.



This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCCID: 2AUMJEV-04-LTE

Support:

1-888-ALOECARE (1-888-256-3772)

support.aloecare.com

support@aloecare.com

Designed by

Aloe Care Health, Inc.

119 West 24th Street, Suite 315

New York, NY 10011, United States

Copyright © 2026



Please Recycle